MCW Chair Weighing Scale





The MCW 300L chair scale provides an affordable and practical solution to address the challenge of weighing patients with limited mobility or who are unable to stand. MCW combines sturdy construction with functionality, minimizing the difficulty of lifting patients and eliminating the need to move heavy health scales. The scale's swivel feet and arm rests add comfort for patients during the weighing process, making it easier for them to get in and out of the chair. Portability is not an issue, as the swivel wheels pivot smoothly, while the scale's built-in battery eliminates the inconvenience of being tethered to an electrical outlet by a power cord. Weight is measured in kilograms or pounds and the MCW 300L calculates body mass index (BMI).

Convenience and accessibility when weighing patients with limited mobility







Features

- Bright LED display for easy viewing in most environments
- · Color-coded keys facilitate quick recognition of the most frequently used buttons
- Splash-proof housing, keypad and display are easy to clean
- · Footrests and armrests pivot for easy patient accessibility
- · Wheels swivel, allowing easy movement
- · BMI (body mass index) indicator
- · Zero-tracking feature ensures display returns to zero reading
- · Rechargeable battery included for operation almost anywhere
- · AC adapter included

Applications

- Weighing
- · BMI (body mass index) indicator

About ADAM:

For more than 45 years, professionals worldwide have relied on Adam Equipment for an extensive selection of dependable and affordable weighing equipment. Adam provides the right balance of speed, performance and value.





Model	MCW 300L
Capacity	660lb / 300kg
Readability	0.1lb / 0.05kg
Repeatability (S.D.)	0.2lb / 0.1kg
Linearity (+/-)	0.2lb / 0.1kg
Weighing Units	kg, lb
Stabilization Time (sec)	3
Interface	RS-232
Body Mass Index [BMI]	Yes
Height Rod	No
Calibration	External calibration
Display	LED with 0.9"/ 22mm -high digits
Power Supply	15VDC 500mA adapter / Internal rechargeable battery
Operating Temperature	5° to 40°C / 41° to 104°F
Overall Dim.	26.6"x34.4"x37.4" / 675x875x950mm (wxdxh)
Net Weight	55.1lb / 25kg

Accessories

Item number	Description
3014011014	RS-232 cable M-F
3074010507	RS-232 to USB adapter (must order corresponding RS-232 cable)
1120014641	AIP impact printer
1120011156	ATP thermal printer
600002028	AdamDU Data Collection Program
700660290	Calibration certificate

Supplied by

www.adamequipment.com